



INTRODUCING INHALED SEDATION

SEDATION IN CRITICALLY ILL PATIENTS

CHAPTER 1

Key role of sedation

Critically ill patients treated with invasive mechanical ventilation very often need to be sedated for their **comfort** and **safety**.¹

COMFORT

- Short- and long-term reduction of anxiety and distress

SAFETY

- Optimising ventilator treatment
- Avoiding self-extubation
- Tolerating procedures
- Reducing autonomic stress

1. Hughes et al., Clin Pharmacol 2012.

Common objectives in sedation recommendations

Sedation recommendations¹⁻³ developed from the Intensive Care medical societies of different countries have common objectives:

Limit long-term physical and psychological impact of ICU stay

Minimize ventilation duration and related complications

Limit side effects using lowest effective dose of sedative (best benefit-risk ratio)

1. Devlin et al., Crit Care Med 2018 (PADIS Clinical Practice Guidelines 2018). 2. Celiz-Rodriguez et al., Med Intensiva (Engl Ed) 2020 (Panamerican Guidelines 2020). 3. DAS-Task Force 2015, Ger Med Sci 2015 (German DAS Guidelines 2015).

Challenges with current sedation strategies

Patient diversity

- Unpredictable PK and PD¹
 - Chronic illness
 - Acute organ dysfunction
 - Advanced age
- Sedation needs change over time
- Development of drug tolerance²

Risk of iatrogenic harm

- Short-term risks of sedation
 - Hypotension, bradycardia²
 - Propofol infusion syndrome³
 - Withdrawal syndrome²
 - Delirium⁴ and hallucinations²
- Long-term risks of sedation
 - Slow and unpredictable awakening, prolonged mechanical ventilation and ICU stay²
 - Cognitive dysfunction⁴
 - Post-traumatic stress disorder²

Finding the right balance

Oversedation^{1,4}
vs
Undersedation⁵

1. Devlin et al., Crit Care Med 2018 (PADIS Clinical Practice Guidelines 2018).
2. Jerath et al., Am J Resp Crit Care Med 2016. 3. Krajcova et al., Crit Care 2015.
4. Pandharipande et al., JAMA 2007. 5. Olsen et al., N Engl J Med 2020.


Consequences of over- and undersedation in mechanically ventilated patients

Undersedation

- Hypertension, tachycardia¹
- Agitated delirium¹
- Patient safety risks^{1,2}
- Ventilator asynchrony¹
- Post-traumatic stress disorder²

Oversedation

- Hypotension¹
- Prolonged ventilator time¹
- Complications to critical care¹
- Delirium³



1. Devlin et al., Crit Care Med 2018 (PADIS Clinical Practice Guidelines 2018). 2. Franks et al., Ann Am Thorac 2021. 3. Pandharipande et al., Anesthesiology 2006.

Proposed strategies to avoid iatrogenic harm

Avoid benzo-diazepines

- Benzodiazepines are associated with the development of delirium¹
- Benzodiazepine infusions lead to prolonged wake-up²
- Guidelines recommend non-benzodiazepine sedation strategies²

Daily wake-up

- Daily interruption of IV sedatives and opioids decreases residual sedation, leading to shorter ventilator time and ICU stay³
- Guidelines recommend SAT together with SBT⁴, unless patient is lightly sedated

Light sedation

- Light IV sedation is as good as daily interruption⁵

SAT=Spontaneous Awakening Test. SBT=Spontaneous Breathing Test.

1. Pandharipande et al., Anesthesiology 2006. 2. Devlin et al., Crit Care Med 2018 (PADIS Clinical Practice Guidelines 2018). 3. Kress et al., N Engl J Med 2000. 4. Girard et al., Lancet 2008. 5. Mehta et al., JAMA 2014.

ABCDEF bundles

Society of Critical Care Medicine: The ABCDEF Bundle elements individually and collectively can help reduce delirium, improve pain management and reduce long-term consequences for adult intensive care unit (ICU) patients.

Assess, prevent and manage pain

Both Spontaneous Awakening Test (SAT)
and Spontaneous Breathing Test (SBT)

Choice of analgesia and sedation


Delirium – assess, prevent, manage

Early mobility and exercise

Family engagement

Effect of ABCDEF bundle adherence

Association between proportional performance of the ABCDEF bundle and patient-related outcomes



ABCDEF bundle performance showed significant and clinically meaningful improvements in outcomes including:

- Survival
- Mechanical ventilation use
- Coma and delirium
- Restraint use
- ICU readmissions
- Post-ICU discharge disposition

The “No/Light sedation” approach

- Motivated by justified concerns regarding iatrogenic harm caused by oversedation.
- Sedation depth has been found to predict increased risk of death, delirium, and delayed time to extubation.¹
- If a patient is awake, comfortable and safe, don't start any sedation.
- Early trials with small sample sizes^{2,3} suggested the feasibility and safety of the no/light sedation approach and hinted some benefits compared to usual care, leading to larger RCTs being performed.

1. Shehabi et al., Crit Care Med 2018. 2. Strøm et al., Lancet 2010. 3. Shehabi et al., Crit Care Med 2013.

Potential challenges with the “No/Light sedation” approach

- Recent publications¹⁻³ highlight the challenge in using a one-size-fits-all approach.
- No/light sedation not suitable for, and/or difficult to achieve in, a significant proportion of patients.^{1,2}
- Three large RCTs¹⁻³ on no/light sedation found no improvements in terms of:
 - Mortality at 90 days
 - Ventilator-, coma/delirium- or ICU-free days
- Concerns regarding patient comfort and safety with no/light sedation approach:
 - Higher risk of self-removal of ET-tube, GI-tube or arterial line²
 - Increased nurse-to-patient ratio needed to manage non-sedated patients in distress²
 - Increased use of restraints on non-sedated patients³

1. Shehabi et al., N Engl J Med 2019. 2. Olsen et al., N Engl J Med 2020 (Electronic supplementary material). 3. Hughes et al., N Engl J Med 2021.

Studies on the “No/Light sedation” approach

Recent publications highlight the challenge in using a one-size-fits-all approach

The NEW ENGLAND JOURNAL of MEDICINE

1

ORIGINAL ARTICLE

Early Sedation with Dexmedetomidine in Critically Ill Patients

Y. Shehabi, B.D. Howe, R. Bellomo, Y.M. Arabi, M. Bailey, F.E. Bass, S. Bin Kadiman, C.J. McArthur, L. Murray, M.C. Reade, I.M. Seppelt, J. Takala, M.P. Wise, and S.A. Webb, for the ANZICS Clinical Trials Group and the SPICE III Investigators*

- No difference in 90-day mortality between dexmedetomidine or usual care
- ~50% of patients *needed* deeper sedation (RASS -3 or deeper) during the first 2 days
- >25% patients needed deeper sedation during the first week

The NEW ENGLAND JOURNAL of MEDICINE

2

ORIGINAL ARTICLE

Dexmedetomidine or Propofol for Sedation in Mechanically Ventilated Adults with Sepsis

C.G. Hughes, P.T. Mailloux, J.W. Devlin, J.T. Swan, R.D. Sanders, A. Anzueto, J.C. Jackson, A.S. Hoskins, B.T. Pun, O.M. Orun, R. Raman, J.L. Stollings, A.L. Kiehl, M.S. Duprey, L.N. Bui, H.R. O'Neal, Jr., A. Snyder, M.A. Gropper, K.K. Guntupalli, G.J. Stashenko, M.B. Patel, N.E. Brummel, T.D. Girard, R.S. Dittus, G.R. Bernard, E.W. Ely, and P.P. Pandharipande, for the MENDS2 Study Investigators*

- No difference in outcomes in patients who received dexmedetomidine compared with propofol
- Despite targeting light sedation a large proportion of patients needed deeper sedation during the first week

The NEW ENGLAND JOURNAL of MEDICINE

3

ORIGINAL ARTICLE

Nonsedation or Light Sedation in Critically Ill, Mechanically Ventilated Patients


Hanne T. Olsen, M.D., Helene K. Nedergaard, M.D., Ph.D., Thomas Strøm, M.D., Ph.D., Jakob Oxlund, M.D., Karl-Andre Wian, M.D., Lars M. Ytrebø, M.D., Ph.D., Bjørn A. Kroken, M.D., Michelle Chew, M.D., Ph.D., Serkan Korkmaz, Jørgen T. Lauridsen, M.Sc., and Palle Toft, M.D., D.M.Sc.

- No difference in 90-day mortality between patients on no sedation compared to light sedation
- Higher risk of self-removal of ET-tube, GI-tube or arterial line in non-sedated patients

1. Shehabi et al., N Engl J Med 2019. 2. Hughes et al., N Engl J Med 2021.
3. Olsen et al., N Engl J Med 2020.

No/Light sedation is not deemed suitable for all patients

Patients with a Clinical Indication for Deep Sedation



No. at Risk

Dexmedetomidine	1952	1915	1775	1551	1351	1151	991	849	747	645	583	515	453	407
Usual care	1963	1928	1798	1610	1384	1201	1045	921	798	698	613	550	496	463

Early Sedation with Dexmedetomidine in Critically Ill Patients

Y. Shehabi, B.D. Howe, R. Bellomo, Y.M. Arabi, M. Bailey, F.E. Bass, S. Bin Kadiman, C.J. McArthur, L. Murray, M.C. Reade, I.M. Seppelt, J. Takala, M.P. Wise, and S.A. Webb, for the ANZICS Clinical Trials Group and the SPICE III Investigators*


- Target sedation level: **RASS -2 to +1**
- Patients in need for deeper sedation (RASS -3 or deeper)
 - First 2 days: ~50%
 - First week: >25%
- 75% of patients on dexmedetomidine needed supplemental sedatives the first two days to achieve the sedation level needed

No/Light sedation might be hard to achieve



Dexmedetomidine or Propofol for Sedation in Mechanically Ventilated Adults with Sepsis

C.G. Hughes, P.T. Mailloux, J.W. Devlin, J.T. Swan, R.D. Sanders, A. Anzueto, J.C. Jackson, A.S. Hoskins, B.T. Pun, O.M. Orun, R. Raman, J.L. Stollings, A.L. Kiehl, M.S. Duprey, L.N. Bui, H.R. O'Neal, Jr., A. Snyder, M.A. Gropper, K.K. Guntupalli, G.J. Stashenko, M.B. Patel, N.E. Brummel, T.D. Girard, R.S. Dittus, G.R. Bernard, E.W. Ely, and P.P. Pandharipande, for the MENDS2 Study Investigators*

- Despite targeting light sedation (RASS 0 to -2), a large proportion of patients required deeper sedation during the first week



No sedation - more risk than benefit?



- Higher risk of self-removal of ET-tube, GI-tube or arterial line in non-sedated patients.
- Self-extubation and removal of other equipment were common **despite a 1:1 nurse-to-patient ratio** in most participating ICUs.

Proposed medical properties of the ideal sedative agent


The 'ideal sedative agent' should possess the following qualities:

- Sedative, analgesic and anxiolytic properties^{1,3}
- Rapid onset and offset of action^{1,3}
- Predictable dose-response relationship¹
- Few side effects¹
- Organ-independent elimination²
- No active metabolites^{2,3}
- No drug accumulation^{1,3}
- Minimal drug interactions^{1,3}

Inhaled sedation (isoflurane/ sevoflurane) has been proposed as "the ideal sedative agent."⁴⁻⁶

1. Ostermann et al., JAMA 2000. 2. Hughes et al., Clin Pharmacol 2012. 3. https://www.ics.ac.uk/Society/Guidance/PDFs/Analgesia_and_Sedation. 4. Spencer et al., Int Care Med 1992. 5. Kong et al., BMJ 1989. 6. Hendrickx et al., J Clin Monit Comput 2018.

Patient-tailored sedation and analgesia



- **UK Analgesia and Sedation guidelines (2014):** "There is insufficient evidence to recommend one regimen over another, and so the agents chosen should be individualized to the patient's requirements, characteristics and the clinical situation."¹
- There is no "one-size-fits-all" solution²

1. https://www.ics.ac.uk/Society/Guidance/PDFs/Analgesia_and_Sedation.
 2. Sackey et al., Anesthesiology 2010.

Current sedative agents


	MIDAZOLAM¹	PROPOFOL²	DEXMEDETOMIDINE³⁻⁵	ISOFLURANE⁶⁻⁹	SEVOFLURANE^{9,10}
Drug class	benzodiazepine	GABA agonist	alpha-2 receptor agonist	fluorinated ether	fluorinated ether
Metabolisation	40-60%	66%	95%	0.2%	2-5%
Degradation metabolite	α-hydroxymidazolam	Glucuronides	Glucuronides	Trifluoroacetic acid, fluorides	Fluorides, hexafluoroisopropanol
Active metabolite	Yes	No	No	No	No
Elimination	60-80% renal	88% renal	90% renal	Pulmonary exhalation	
Duration of action T _{1/2} context-sensitive	1.5-3.0 h (single shot)	1.5-12.4 h	1.9-2.5 h	< 5 min	< 5 min
80% decrement time	3-12 h (up to several days)	3.5 h-3 d	3.7 h	30-35 min	< 8 min

1. Nordt et al., J Emerg Med 1997. 2. Trapani et al., Curr Med Chem 2000. 3. Weerink et al., Clin Pharmacokin 2017. 4. Iirola et al., Crit Care 2011. 5. Venn et al., Br J Anaesth 2002. 6. Mazze et al., Anesthesiology 1974. 7. Holaday et al., Anesthesiology 1975. 8. Kharasch et al., Anesthesiology 1999. 9. Bailey et al., Anesth Analg 1997. 10. Behne et al., Clin Pharmacokin 1999.

Pharmacokinetics of intravenous vs inhaled sedatives

Elimination of intravenous sedatives relies on adequate hepatic and renal function¹

Hepatic and renal function are often impaired in ICU patients², leading to slow elimination¹, risk of drug accumulation and oversedation³⁻⁸




Elimination of inhaled sedatives is independent of hepatic and renal function^{9,10}

Inhaled sedation works even in severe ARDS on ECMO^{11,12}

1. Jerath et al., Anesth Analg 2017. 2. Vincent et al., Crit Care Med 2006. 3. Shelly et al., Eur J Anaesthesiol 1991. 4. Devlin et al., Crit Care Med 2018 (PADIS Clinical Practice Guidelines 2018). 5. Barr et al., Anesthesiology 2001. 6. Jerath et al., 2016. 7. Shafer et al., Crit Care Med 1998. 8. Sackey et al., Crit Care Med 2004. 9. Eger, Br J Anaesth 1984. 10. Lu et al., Pharmacology 2008. 11. Meiser et al., Respir Care 2018. 12. Rand et al., J Artif Organs 2018.


History of Inhaled Sedation



*Currently called Sedaconda® ACD (Anaesthetic Conserving Device)

1. Kong et al., BMJ 1989. 2. Meiser et al., Lancet Resp Med 2021.

The Sedaconda study



- The largest prospective trial to date on inhaled sedation
- A phase 3, randomised, controlled, open-label, multicentre, parallel-group, non-inferiority trial designed to meet the European regulatory requirements for approval
- The study objective was to demonstrate the efficacy and safety of isoflurane for sedation in invasively ventilated intensive care unit (ICU) patients using the Sedaconda ACD delivery device
- The study demonstrated non-inferior sedation efficacy of isoflurane compared to propofol


Isoflurane for Inhaled Sedation in the ICU

- Sedaconda® (isoflurane) is currently the only drug approved for the indication of inhaled sedation of adults in the ICU.¹
- Metabolised at a very low rate (0.2%)²
 - Elimination independent of liver or renal function³
- Tissue accumulation in muscle and fat during sedation is clinically insignificant⁴
- Elimination is via the lungs (>99%)^{3,5,6}
 - Wake-up time within one hour after days of sedation^{4,7,8}

Pharmacodynamics - How do inhaled sedatives work?


Multiple mechanisms of action¹⁻⁴

- Potentiation of inhibitory receptors
- Inhibition of excitatory receptors
- Direct cell membrane effects



1. Campagna et al., N Engl J Med 2003. 2. Herold et al., J Gen Physiol 2014.
3. Wu et al., Anesthesiology 2004. 4. Herold et al., Eur Biophys J 2017.

Pharmacodynamics of isoflurane, besides sedation



1. Crystal, J Anesth Hist 2017. 2. Turner et al., Respir Care 2012. 3. Basil et al., Anesthesiology 1999.

EFFICACY OF INHALED SEDATION

CHAPTER 2

Inhaled sedatives – Fast onset of sedation

Sedation efficacy



**Administered
via the airways
and lungs**



**Rapid
alveolar
uptake**



**Delivered via
the bloodstream
to the brain**




**Fast
onset of
sedation**

- No significant concerns of drug tolerance
- Dose requirements do normally not increase over time

Sedation efficacy isoflurane vs propofol

The Sedaconda study:

Comparable time spent in the target RASS range without rescue sedation




Proportion of time within the sedation target range:

- >90% in both groups
- Isoflurane non-inferior to propofol

RASS=Richmond Agitation-Sedation Scale

Lower opioid requirements with isoflurane


Opioid reduction



MED=Morphine Equivalent Dose. BPS=Behavioural Pain Scale

- In the Sedaconda study, opioid dose intensity was **29% lower** in the isoflurane group compared to the propofol group.
- No signs of clinical pain issues in relation to opioid dose reduction, as indicated by BPS scores.

Opioid reduction in Covid-19 ARDS after switching from IV sedation to isoflurane




- Switch from midazolam to isoflurane in 11 mechanically ventilated Covid-19 patients.
- Sufentanyl consumption decreased significantly while the same analgo-sedation goal was reached.

Opioid reduction with isoflurane¹

- Inhaled anaesthetics have antinociceptive effects on the spinal cord^{1,2} which may explain the reduced opioid needs with isoflurane.
- Lower opioid doses may lead to:
 - Gut mobility better preserved³
 - Spontaneous breathing preserved¹
 - Risk of delirium decreased^{4,5}

Spontaneous breathing more common with isoflurane



Number of patients	
Isoflurane	149 148 146 142 117 93 87 93 95 92 80 33
Propofol	151 151 150 146 129 91 89 93 93 89 79 44

Estimated rate of spontaneous breathing in the Sedaconda study:

Day 1

- Isoflurane 50.3%
 - Propofol 37.0%
- } p=0.013

Day 2

- The difference favored isoflurane but did not reach statistical significance

Spontaneous breathing preserves muscle function

- Spontaneous breathing prevents respiratory muscle fiber dysfunction^{1,2}
- Preserved respiratory muscle function shortens ICU length of stay³
- Respiratory muscle weakness after mechanical ventilation is associated with one-year mortality⁴
- Augmented spontaneous breathing may recruit lung tissue in dependent areas⁵


“Improvements in pulmonary gas exchange, systemic blood flow and oxygen supply to the tissue which have been observed when spontaneous breathing has been maintained during mechanical ventilation are reflected in the clinical improvement in the patient’s condition.”³

Isoflurane facilitated opioid reduction and spontaneous breathing in ARDS patients

Opioid reduction

Spontaneous breathing

Opioid consumption before and during isoflurane sedation compared with propofol/midazolam



Data are shown as mean ± SD.


*P compared to baseline, #P compared to propofol/midazolam.

	Propofol/ Midazolam n = 19	Isoflurane n = 19	P
SAPS II, points			
Before	43.2 ± 15.2	40.2 ± 9.6	.47
6 h	41.4 ± 14.9	39.2 ± 9.8	.61
24 h	42.6 ± 13.8	35.7 ± 10.2	.10
Remifentanyl, µg/kg/min			
Before	n = 14	n = 16	
6 h	0.22 ± 0.09	0.19 ± 0.10	.39
24 h	0.23 ± 0.10	0.10 ± 0.04*	.007
24 h	0.25 ± 0.09	0.09 ± 0.04*	< .001
Sufentanil, µg/kg/h			
Before	n = 5	n = 3	
6 h	0.68 ± 0.59	0.46 ± 0.66	.64
24 h	0.68 ± 0.58	0.29 ± 0.45	.20
24 h	0.52 ± 0.55	0.29 ± 0.45	.38
Spontaneous breathing			
Before	3 (16)	2 (11)	.64
6 h	3 (16)	12 (63)	.003
24 h	3 (16)	17 (90)	< .001

Shorter and more predictable wake-up time with isoflurane

Short time to wake-up

Time to wake-up during SAT




SAT = spontaneous awakening trial.

The Sedaconda study showed:

- Faster wake-up with isoflurane than propofol
 - Day 1
 - No statistically significant difference
 - Day 2
 - Isoflurane: 20 min (IQR 10-30 min)
 - Propofol: 30 min (IQR 11-120 min)
- Lower inter-individual variability with isoflurane (more predictable)

Shorter ventilator time and ICU stay with isoflurane in cardiac arrest patients

Shorter stay in ICU



Isoflurane sedation was associated with:

- Shorter ventilator time (p=0.003)
- Shorter ICU stay (p=0.006)

Fast cognitive recovery

Fast cognitive recovery

1

Isoflurane compared with midazolam for sedation in the intensive care unit

Kin Leong Kong, Sheila M Willatts, Cedric Prys-Roberts

- The median times from stopping sedation to patients writing their home address were significantly shorter for patients sedated with isoflurane compared with those sedated with midazolam¹

2

Prolonged isoflurane sedation of intensive care unit patients with the Anesthetic Conserving Device

Peter V. Sackey, MD; Claes-Roland Martling, MD, PhD; Fredrik Granath, PhD; Peter J. Radell, MD, PhD

- Patients sedated with isoflurane emerged from sedation faster than after midazolam²
- Time to follow verbal command²
 - Isoflurane: 10 ± 8 min
 - Midazolam: 110 ± 130 min

3

Short- and long-term follow-up of intensive care unit patients after sedation with isoflurane and midazolam—A pilot study*


Peter V. Sackey, MD, PhD; Claes-Roland Martling, MD, PhD; Christine Carlswård, MD; Örjan Sundin, PhD; Peter J. Radell, MD, PhD

- Sedation of ICU patients with isoflurane may result in fewer delusional memories or hallucinations from the ICU compared with more commonly used intravenous sedation³

1. Kong et al., BMJ 1989. 2. Sackey et al., Crit Care Med 2004.
3. Sackey et al., Crit Care Med 2008.


Responsiveness after sedation

Fast cognitive recovery



Fast cognitive recovery after 96 hours of sedation

Time to follow verbal command from termination of sedation



- Patients sedated with isoflurane emerged from sedation faster than with midazolam
- Time to follow verbal command
 - Isoflurane: 10 ± 8 min
 - Midazolam: 110 ± 130 min

} p=0.003

Isoflurane sedation in ICU patients has been associated with a low incidence of hallucinations or delusions

Fast cognitive recovery



Sackey et al., Crit Care Med 2008.

M/EN/PH/220001 Nov 2021

Rapid and predictable emergence from sedation is a clinically valuable feature of a sedative in the ICU

Facilitates **planning** of extubation¹ and **aftercare**, including mobilization^{2,3}

Facilitates reliable **neurological evaluation**^{4,5} and reduces the need for diagnostic tests⁶

Enables **engagement** with family and staff²⁻⁵



1. Meiser et al., Lancet Resp Med 2021. 2. Schweickert et al., Lancet 2009.
3. Green et al., J Multidiscip Healthc 2016. 4. Kong et al., BMJ 1989. 5. Sackey et al., Crit Care Med 2004. 6. Kress et al., N Engl J Med 2000.

EXPLORATORY STUDIES & FINDINGS

CHAPTER 3


Preclinical observations in ARDS animal models (inhaled vs IV sedation)

- Reduced levels of inflammatory biomarkers¹⁻⁴
- Reduced levels of markers for epithelial and endothelial injury¹⁻⁴
- Improved arterial oxygenation observed as ratio of oxygen tension to inspired oxygen fraction ($\text{PaO}_2/\text{FiO}_2$)¹⁻⁴

ARDS = Acute Respiratory Distress Syndrome

Clinical observations in ARDS patients

Evolution of PaO₂/FiO₂ ratio (mmHg) during the 120-hour observation period (intervention during the first 48 hours)



- On Day 2, PaO₂/FiO₂ ratio was higher in the sevoflurane group than in the midazolam group
- On Day 2, both plasma and alveolar sRAGE levels (an epithelial injury marker) were lower in the sevoflurane group than in the midazolam group

sRAGE=soluble form of the receptor for advanced glycation end-products

ARDS, cytokines and other indicators - Effects of inhaled anaesthetics

EFFECTS OF INHALED ANAESTHETICS ON ARDS		General findings in ARDS (Sweeney 2016) ¹	Findings in Corona virus infection and ARDS (Huang 2020) ²	Effects of Inhaled Anaesthetics (isoflurane and/or sevoflurane) vs Intravenous Anaesthetics (controls)				
				Voigtsberger 2009 ³ (rat model)	Ferrando 2013 ⁴ (pig model)	Strosing 2016 ⁵ (mouse model)	Kellner 2017 ⁶ (rat model)	Jabaudon 2017 ⁷ (human)
Pro- inflammatory cytokines	TNF-alpha	↑	↑	↓	↓	NA	NA	↓
	IL-1 beta	↑	↑	↓	↓	↓	NA	NS
	IL-6	↑	NA	↓	↓	NA	↓	↓
	IL-8	↑	↑	NA	↓	NA	NA	↓
	IL-10	↑	↑	NA	NA	NA	NS	NA
	MCP-I	↑	↑	↓	NA	NA	NS	NA
	CINC-I	↑	NA	↓	NA	NA	↓	NA
Other markers of lung injury	White blood cell count in alveoli	↑	NA	↓	↓	↓	↓	NA
	Lung fluid permeability/edema	↑	NA	↓	↓	NA	↓	NA
	Alveolar histological disruption	↑	NA	Maintained	NA	Maintained	NA	NA
	S-RAGE	↑	NA	NA	NA	NA	NA	↓
	Oxygenation	↓	↓	↑	↑	↑	↑	↑

ARDS: Acute Respiratory Distress Syndrome; NA: Not assessed; NS: No significant difference

1. Sweeney et al., Lancet 2016. 2. Huang et al., Lancet 2020. 3. Voigtsberger et al. Anesthesiol 2009. 4. Ferrando et al., Eur J Anaesthesiol 2013. 5. Strosing et al., Anesth Analg. 2016. 6. Kellner et al., Anesth Analg. 2017. 7. Jabaudon et al., Am J Resp Crit Care Med 2016. 8. Sackey et al., Crit Care Med. 2004. 9. Mesnil et al., Intensive Care Med. 2011.

Ongoing clinical studies (ARDS)

The SESAR study¹

Patient population:

700 intensive care patients with ARDS

Study treatment:

Inhaled sedation with sevoflurane for up to seven days vs intravenous propofol

Primary objective:

Efficacy of inhaled sedation measured as ventilator-free days at day 28, considering death as a competing event

The ISCA study²

Patient population:

400 COVID-19-related ARDS patients

Study treatment:

Inhaled sedation (isoflurane or sevoflurane) versus intravenous sedation

Primary objective:

Efficacy of inhaled sedation measured as ventilator-free days at day 28, considering death as a competing event

The SAVE-ICU study³

Patient population:

a) Patients with proven or suspected COVID-19, or b) COVID-19 negative patients who have a PaO₂/FiO₂ ratio ≤300 measured with arterial blood gas at least once during the 12 hours prior to enrollment

Study treatment:

Inhaled sedation (isoflurane or sevoflurane) vs standard care (any iv sedation)

Primary objective:

Efficacy of inhaled sedation measured as 1) hospital mortality within 2 years, 2) ventilator-free days within 30 days, 3) ICU-free days within 30 days, 4) participant Quality of Life at 3 and 12 months after discharge.

ARDS = Acute Respiratory Distress Syndrome

ClinTrials.gov Identifiers: 1. NCT04235608. 2. NCT04383730. 3. NCT04415060.

Reduced bronchospasm

Pulmonary protection

Intensive Care Med (2006) 32:927–933
DOI 10.1007/s00134-006-0163-0

PEDIATRIC BRIEF REPORT

Venkat Shankar
Kevin B. Churchwell
Jayant K. Deshpande

Isoflurane therapy for severe refractory status asthmaticus in children

Isoflurane resulted in:

- An immediate clinical improvement in all 11 children
- An improvement in arterial pH
- A reduction in partial pressure of arterial carbon dioxide (PaCO_2) in all 11 patients

Rapid weaning from mechanical ventilation occurred in 10 patients.


Isoflurane for Life-Threatening Bronchospasm: A 15-Year Single-Center Experience

David A Turner MD, David Heitz RRT, Mehrengise K Cooper FRCPH,
P Brian Smith MD MPH MHS, John H Arnold MD, and Scot T Bateman MD

Isoflurane led to improvement in pH and pCO_2 in patients with life-threatening bronchospasm:

- Between 4 to 24 hours, there was a statistically significant decrease in PIP ($p=0.006$)
- FIO_2 decreased within 4 hours of initiation of isoflurane ($p=0.001$)
- FIO_2 decreased from 4 to 24 hours ($p =0.02$)

Ongoing clinical studies on Inhaled Sedation in intensive care



SAFETY OF INHALED SEDATION

CHAPTER 4

Safety of Inhaled Sedation with isoflurane

Adverse reactions to isoflurane are generally dose-dependent extensions of pharmacologic effects¹:

- hypotension
- respiratory depression



1. Sedaconda® SmPC – August 2021. 2. Meiser et al., Lancet Resp Med 2021.

Important safety information / Precautions

Malignant Hyperthermia (MH)

- Sedation with isoflurane is contraindicated in patients with known or suspected genetic susceptibility to MH.¹
- MH is a rare genetic disorder (incidence 1/10.000–250.000²) where isoflurane sedation may trigger a skeletal muscle hypermetabolic state.

Intracranial Pressure (ICP)

- During sedation with isoflurane, ICP may increase slightly.¹
- Caution should be taken when administering isoflurane to patients with increased ICP, and ICP must be monitored in such patients.¹

1. Sedaconda® SmPC – August 2021. 2. Rosenberg et al., Orphanet J Rare Dis 2015.

Malignant Hyperthermia (MH)


- In susceptible individuals, isoflurane sedation may trigger a skeletal muscle hypermetabolic state, leading to high oxygen demand and the clinical syndrome known as malignant hyperthermia (MH).
- Known or suspected genetic susceptibility to MH is a contraindication to inhaled sedation.¹
- Low incidence of MH reactions (1:10,000 - 1:250,000).²
- The syndrome includes non-specific features such as:
 - muscle rigidity, tachycardia, tachypnoea, cyanosis, arrhythmias, and unstable blood pressures
- Treatment of MH includes:
 - discontinuation of triggering agent (isoflurane)
 - IV administration of antidote dantrolene sodium
 - supportive therapy

Malignant Hyperthermia and PRIS at a glance

	Malignant Hyperthermia (MH)¹⁻⁴	Propofol-infusion syndrome (PRIS)⁵⁻⁸
Triggers	Volatile anaesthetics Succinylcholine	Propofol
Dose-dependent	No	Yes (>4 mg/kg/d; >48h)
Incidence	1:10.000-1:250.000	1,1%
Therapy	Antidote dantrolene Withdrawal of agent, cooling	No antidote Symptomatic
Prevention	Medical history of family MH No use in muscle disease	<4 mg/kg/h, <7 d No use in children <16 y
Mortality	<5%	51%

1. Larach et al., Anesthesiology 1994. 2. Glahn et al., Brit J Anaesth 2010. 3. Rosenberg et al., Orphanet J Rare Dis 2015. 4. Bandschapp et al., Swiss Med Wkly 2012. 5. Roberts et al., Critical Care 2009. 6. Krajčová et al., Critical Care 2015. 7. Mirrakhimov et al., Crit Care Res Pract 2015. 8. Eziefule et al., Am J Perinatol Rep 2016.

Low risk of room pollution and exposure of anaesthetic gases to the ICU staff



1. Sackey et al., Crit Care 2004. 2. Gonzalez-Rodriguez et al., Rev Esp Anesthesiol Reanim 2014. 3. Pickworth et al., Can J Anaesth 2013. 4. Contreras et al., ESICM Congress abstract 2021.

ELIGIBLE PATIENTS FOR INHALED SEDATION

CHAPTER 6

Advantages of Inhaled Sedation

- Reduces the need for multiple intravenous sedatives¹⁻³
- Reduces opioid doses, potentially reducing opioid side effects¹⁻⁴
- Facilitates spontaneous breathing^{3,4}
- Minimises residual sedation and leads to rapid and predictable wake-up⁴ and return of cognitive function^{5,6}

When can inhaled sedation be a good choice for mechanically ventilated ICU patients

